




Holly End House  
Water Lane  
Bovingdon  
Herts  
HP3 0NA  
T: 00 44 (0) 1442 832764  
Fx: 00 44 (0) 1442 832550  
Em: [info@inagas.com](mailto:info@inagas.com)  
[www.inagas.com](http://www.inagas.com)

Producers of insulated glass who have passed the EN 1279 parts 2 & 3 for gas filled units are required to test gas concentration levels on a regular basis to maintain their certification and ensure continued consistency of quality.



The latest Testoxy instrument is easy to use and powered by maintenance free batteries which can be recharged using the provided charger. The Testoxy can be used with all I.G gasses such as Argon, Krypton, Xenon and pre-mixed gasses.

The testing is done by simply inserting the probe into the sealed unit and pressing the start button. The LED display starts by showing the Inagas logo which indicates the warm up and calibration sequence has started. The screen will then display 00 and is ready to take a gas sample, the result is then shown as a percentage of gas in the IG unit.

#### Specifications

Dims: 140 x 85 x 45mm

Weight: 0.4kg

Power: Alkaline AA Cells

NiMH recharge pack

#### Supplied with the following

Sampling hose

Re-Charging unit

Piercing tool

Instructions on use

Calibration certificate